

Dear,

I hope you are keeping well in these strange times. Even if we are not allowed to visit you, we believe it is of outmost importance to let you know that we are here to support you.

With this Newsletter, we would like to give you an update about the MYNX CONTROL™ Vascular Closure Device (VCD) – our next-generation of delivery. A tension indicator provides visual confirmation of device positioning for effective deployment, and much more. Below you will find some more information that I hope you will find interesting. Please don't hesitate to contact me for further information. Take care of yourself!

With very best regards,

Kia Lovberg

Field Product Manager Nordics, Cordis

+46 703 12 22 38 / katarina.lovberg@cardinalhealth.com

MYNX CONTROL™ Vascular Closure Device CLOSES WITH SECURITY. LEAVES WITHOUT A TRACE.

MYNX CONTROL™
VASCULAR CLOSURE DEVICE

Closes with Security. Leaves without a Trace.

SECURE CLOSURE **SAFETY AND PATIENT COMFORT** **EASE OF USE**

Next-generation delivery system

- Deploy and retract: 2-button operation - simultaneous procedure
- Sheath catch supports consistent deployment
- Sealant is resorbed by the body within 30 days

Tension indicator provides visual confirmation of device positioning for effective deployment

- Visual confirmation on handle
- Correct balloons-to-wall apposition, for effective deployment

Deploys through standard procedural sheath

- No sheath exchange or up-sizing required

Further information on MYNX CONTROL™ VCD [click here](#)

CE

Note: MYNX CONTROL™ VCD is incompatible with Medtronic Input® Introducer (11cm) sheaths, Cook Check-Flo® Performer® Introducer sheaths, and procedural sheaths longer than 12cm in effective length.

For Healthcare Professionals Only.

Important information: Prior to use, refer to the instructions for use supplied with this device for indications, contraindications, side effects, suggested procedure, warnings and precautions. As part of its continuous product development policy, Cordis reserves the right to change product specifications without prior notification. Please contact your Cordis representative for additional product availability information.

Please contact your MYNX CONTROL™ representative for additional product availability information. CORDIS, CORDIS Logo and MYNX CONTROL are trademarks of Cardinal Health and may be registered in the US and/or in other countries. © 2021 Cardinal Health. All Rights Reserved. 100581571 04/2021

Please Note: MYNX CONTROL™ VCD is incompatible with Medtronic Input® Introducer (11cm) sheaths, Cook Check-Flo® Performer® Introducer sheaths, and procedural sheaths longer than 12cm in effective length.

For Healthcare Professionals Only.

Important information: Prior to use, refer to the instructions for use supplied with this device for indications, contraindications, side effects, suggested procedure, warnings and precautions. As part of its continuous product development policy, Cordis reserves the right to change product specifications without prior notification. Please contact your Cordis representative for additional product availability information.

If you do not wish to receive any further emails from Cordis about our educational events and products/services, reply to this email with the word 'unsubscribe' in the subject line. More information about how Cordis processes personal data can be found at https://www.cardinalhealth.co.uk/en_gb/search.html?q=privacy%20policy.

CORDIS, CORDIS Logo and MYNX CONTROL are trademarks of Cardinal Health and may be registered in the US and/or in other countries. © 2021 Cardinal Health. All Rights Reserved. 100581571 04/2021