



# SABER™ .035

PTA DILATATION CATHETER

DESIGNED FOR  
**EFFICIENCY**

**Cordis**™  
A Cardinal Health company

# EXCELLENT NAVIGATION THROUGH TORTUOUS ANATOMY AND DIFFICULT LESIONS

## KEY FEATURES FOR ENHANCED PERFORMANCE



Smooth tapered tip enhances crossability



Dual layer hydrophilic coating improves trackability




Dual distal markers on  $\geq 100$  mm balloons enhance visualization under fluoroscopy



Proprietary balloon design offers excellent rewrap and facilitates smooth retraction through the sheath

## BALLOON CROSSABILITY RESULTS<sup>1</sup>



6 x 200 mm Superficial Femoral Artery (SFA) Crossing Model

Note: The data provided is based on bench testing performed in 2020 by Cordis Corporation

# OPTIMAL DESIGN FOR STREAMLINED CARE IN MORE PATIENTS

## EXTENSIVE SIZE MATRIX EXPANDS TREATMENT APPLICATIONS

**SABER™ .035 catheter size matrix chart\***  
Usable Catheter Lengths: 40, 80, 135 cm

BALLOON DIAMETER (mm)	BALLOON LENGTH (mm)									
	20	40	60	80	100	120	150	200	250	300
3	24	24	24	22	22	22	22	22	20	20
4	22	22	22	22	22	22	22	22	20	20
5	22	22	22	22	22	20	20	20	18	18
6	18	18	17	17	17	17	17	17	17	17
7	18	18	16	16	16	14	14	14	14	14
8	16	16	16	16	16					
9	13	13	13	13	13					
10	14	14	14	14	14					
12	12	12								

5F
  6F
  7F

Values in table represent RBP

\*Not all balloon sizes are available in all shaft lengths. Please consult the product catalog for additional information.

## PORTFOLIO THAT ENHANCES RESULTS

◀ 300 mm ▶

### BALLOON LENGTHS OF UP TO 300 mm

enable treatment of long lesions with fewer inflations and as a result may potentially reduce procedural and fluroscopy time

### LOW PROFILE 5F SHEATH COMPATIBILITY

for most sizes may reduce the potential for access site complications


### HIGH RATED BURST PRESSURE

of up to 24 atm for effective dilation of tough lesions

### COMPREHENSIVE PRODUCT PORTFOLIO

reduces the need for multiple vendors

# SABER™ .035 PTA DILATATION CATHETER OPTIMIZES OUTCOMES WITH:



**SUPERIOR  
CROSSABILITY<sup>1</sup>**



**EXCELLENT  
DELIVERABILITY**



**BROAD SIZE  
MATRIX**



Contact your Cordis representative to evaluate SABER™ .035 catheters or for more information on the full line of lower-extremity solutions, visit [www.cardinalhealth.co.uk](http://www.cardinalhealth.co.uk)

**REFERENCES** 1. Cordis internal data on file

For Healthcare Professionals Only.

**Important Information:** Prior to use, refer to the instructions for use supplied with this device for indications, contraindications, side effects, suggested procedure, warnings and precautions. As part of its continuous product development policy, Cordis reserves the right to change product specifications without prior notification.

Please contact your Cordis representative for additional product availability information.

CORDIS, Cordis LOGO and SABER are trademarks of Cardinal Health and may be registered in the US and/or in other countries. All other trademarks are the property of their respective owners.

© 2020 Cardinal Health. All Rights Reserved. [100568106] 09/2020

**Cordis**<sup>™</sup>  
A Cardinal Health company