


# SABER™.035

PTA DILATATION CATHETER

## PRODUCT DESCRIPTION

**TYPE:** Low profile, high pressure, over-the-wire PTA balloon catheter designed for use in the iliac, femoral, ilio-femoral, popliteal, infra popliteal and renal arteries and for the treatment of obstructive lesions of native or synthetic arteriovenous dialysis fistulae. The device is also indicated for post-dilatation of balloon expandable and self-expanding stents in the peripheral vasculature.

<b>RBP RANGE</b>	12-24 atm	<b>SHEATH FIT RANGE</b>	5F-7F
<b>NOMINAL PRESSURE</b>	8-10 atm	<b>BALLOON MATERIAL</b>	Nylon (Grilamid)
<b>BALLOON LENGTHS</b>	20-300 mm	<b>MARKERS</b>	Proximal and distal marker bands (dual distal marker bands for balloon lengths ≥100 mm)
<b>BALLOON DIAMETERS</b>	3-12 mm	<b>GUIDEWIRE</b>	.035"
<b>SHAFT LENGTHS</b>	40 cm, 80 cm, 135 cm	<b>COATING</b>	Dual layer hydrophilic coating
<b>SHAFT MATERIAL</b>	Pebax		



## ORDERING INFORMATION

SABER™ .035 PTA DILATATION CATHETER						
BALLOON		RBP	SHEATH	SHAFT LENGTH	SHAFT LENGTH	SHAFT LENGTH
DIAMETER (mm)	LENGTH (mm)					
3	20	24	5	48003002P	48003002N	48003002R
3	40	24	5	48003004P	48003004N	48003004R
3	60	24	5	48003006P	48003006N	48003006R
3	80	22	5	48003008P	48003008N	48003008R
3	100	22	5	48003010P	48003010N	48003010R
3	120	22	5		48003012N	48003012R
3	150	22	5		48003015N	48003015R
3	200	22	5		48003020N	48003020R
3	250	20	5		48003025N	48003025R
3	300	20	5		48003030N	48003030R

SABER™ .035 PTA DILATATION CATHETER						
BALLOON		RBP	SHEATH	SHAFT LENGTH	SHAFT LENGTH	SHAFT LENGTH
DIAMETER (mm)	LENGTH (mm)					
4	20	22	5	48004002P	48004002N	48004002R
4	40	22	5	48004004P	48004004N	48004004R
4	60	22	5	48004006P	48004006N	48004006R
4	80	22	5	48004008P	48004008N	48004008R
4	100	22	5	48004010P	48004010N	48004010R
4	120	22	5		48004012N	48004012R
4	150	22	5		48004015N	48004015R
4	200	22	5		48004020N	48004020R
4	250	20	5		48004025N	48004025R
4	300	20	5		48004030N	48004030R

# SABER™ .035

PTA DILATATION CATHETER

## ORDERING INFORMATION

SABER™ .035 PTA DILATATION CATHETER						
BALLOON						
DIAMETER (mm)	LENGTH (mm)	RBP	SHEATH	SHAFT LENGTH 40 cm	SHAFT LENGTH 80 cm	SHAFT LENGTH 135 cm
5	20	22	5	48005002P	48005002N	48005002R
5	40	22	5	48005004P	48005004N	48005004R
5	60	22	5	48005006P	48005006N	48005006R
5	80	22	5	48005008P	48005008N	48005008R
5	100	22	5	48005010P	48005010N	48005010R
5	120	20	5	48005012P	48005012N	48005012R
5	150	20	5		48005015N	48005015R
5	200	20	5		48005020N	48005020R
5	250	18	5		48005025N	48005025R
5	300	18	5		48005030N	48005030R
6	20	18	5	48006002P	48006002N	48006002R
6	40	18	5	48006004P	48006004N	48006004R
6	60	17	5	48006006P	48006006N	48006006R
6	80	17	5	48006008P	48006008N	48006008R
6	100	17	5	48006010P	48006010N	48006010R
6	120	17	5	48006012P	48006012N	48006012R
6	150	17	5		48006015N	48006015R
6	200	17	6		48006020N	48006020R
6	250	17	6		48006025N	48006025R
6	300	17	6		48006030N	48006030R
7	20	17	5	48007002P	48007002N	48007002R
7	40	17	5	48007004P	48007004N	48007004R
7	60	16	5	48007006P	48007006N	48007006R
7	80	16	5	48007008P	48007008N	48007008R

SABER™ .035 PTA DILATATION CATHETER						
BALLOON						
DIAMETER (mm)	LENGTH (mm)	RBP	SHEATH	SHAFT LENGTH 40 cm	SHAFT LENGTH 80 cm	SHAFT LENGTH 135 cm
7	100	16	6	48007010P	48007010N	48007010R
7	120	14	6		48007012N	48007012R
7	150	14	6			48007015R
7	200	14	6			48007020R
7	250	14	6			48007025R
7	300	14	6			48007030R
8	20	16	6	48008002P	48008002N	48008002R
8	40	16	6	48008004P	48008004N	48008004R
8	60	16	6	48008006P	48008006N	48008006R
8	80	16	6	48008008P	48008008N	48008008R
8	100	16	6	48008010P	48008010N	48008010R
9	20	13	6	48009002P	48009002N	48009002R
9	40	13	6	48009004P	48009004N	48009004R
9	60	13	6	48009006P	48009006N	48009006R
9	80	13	6	48009008P	48009008N	48009008R
9	100	13	6	48009010P	48009010N	48009010R
10	20	14	7	48010002P	48010002N	48010002R
10	40	14	7	48010004P	48010004N	48010004R
10	60	14	7	48010006P	48010006N	48010006R
10	80	14	7	48010008P	48010008N	48010008R
10	100	14	7	48010010P	48010010N	48010010R
12	20	12	7	48012002P	48012002N	48012002R
12	40	12	7	48012004P	48012004N	48012004R

Contact your Cordis representative to evaluate SABER™ .035 catheters or for more information on the full line of lower-extremity solutions, visit [www.cardinalhealth.co.uk](http://www.cardinalhealth.co.uk)

For Healthcare Professionals Only.

**Important Information:** Prior to use, refer to the instructions for use supplied with this device for indications, contraindications, side effects, suggested procedure, warnings and precautions. As part of its continuous product development policy, Cordis reserves the right to change product specifications without prior notification.

Please contact your Cordis representative for additional product availability information.

CORDIS, Cordis LOGO and SABER are trademarks of Cardinal Health and may be registered in the US and/or in other countries. All other trademarks are the property of their respective owners.

© 2020 Cardinal Health. All Rights Reserved. [100568714] 09/2020

