

Hej,

Jeg håber, du har det godt i denne mærkelige tid. Selvom vi ikke må besøge dig, synes vi, det er meget vigtigt at fortælle dig, at vi stadig er her for dig.

Med dette nyhedsbrev vil vi gerne give dig sidste nyt omkring MYNX CONTROL™ Vascular Closure Device (VCD) (karlukkingsenhed) – vores næste generation af indføringsudstyr. En spændingsindikator sørger for visuel bekræftelse af anordningens positionering med henblik på effektiv placering og meget mere. Nedenfor får du nogle informationer, som jeg håber, du vil finde interessante.

Du er velkommen til at kontakte mig for at få yderligere information. Pas godt på dig selv!

Med venlig hilsen,

Christina Carlsson

Account Manager Sweden South, Cordis

+46 702 20 63 94 / christina.carlsson@cardinalhealth.com

## MYNX CONTROL™ Vascular Closure Device CLOSES WITH SECURITY. LEAVES WITHOUT A TRACE.

**MYNX CONTROL™**  
VASCULAR CLOSURE DEVICE

**Cordis**  
A Cardinal Health company

**CLOSES WITH SECURITY.  
LEAVES WITHOUT A TRACE.**

**SECURE CLOSURE** **SAFETY AND PATIENT COMFORT** **EASE OF USE**

**Next-generation delivery system**  
- Deploy and retract: 2-button operation streamlines procedure  
- Sheath catch supports consistent deployment  
- Sealant is resorbed by the body within 30 days

**Tension indicator provides visual confirmation of device positioning for effective deployment**  
- Visual confirmation on handle  
- Correct balloon-to-wall apposition, for effective deployment

**Deploys through standard procedural sheath\***  
- No sheath exchange or up sizing required

Further information on MYNX CONTROL™ VCD click here.

**Step 1: Position the Balloon**  
Result: Temporary Hemostasis

**Step 2: Deploy the Sealant**  
Result: Delivery

**Step 3: Remove the Device**  
Result: Secure Extravascular Closure

CE 2797

\*New MYNX CONTROL™ VCD is incompatible with Medtronic Input® Introducer (11cm) sheaths, Cook Check-Flo® Performer® Introducer sheaths, and procedural sheaths longer than 12cm in effective length.

© Cordis Corporation 2021. All rights reserved. Information that is not in the instructions for use supplied with this device for indications, contraindications, side effects, suggested procedure, warnings and precautions. Part of our continuous product development policy, Cordis reserves the right to change product specifications without prior notification. Please contact your Cordis representative for additional product availability information. CORDIS, CORDIS Logo and MYNX CONTROL are trademarks of Cordis Health and are registered trademarks of Cordis Health in the United States and other countries. All other trademarks are the property of their respective owners.

© 2021 Cordis Health. All Rights Reserved. 1861111 01/2021

Please Note: MYNX CONTROL™ VCD is incompatible with Medtronic Input® Introducer (11cm) sheaths, Cook Check-Flo® Performer® Introducer sheaths, and procedural sheaths longer than 12cm in effective length.

For Healthcare Professionals Only.

Important information: Prior to use, refer to the instructions for use supplied with this device for indications, contraindications, side effects, suggested procedure, warnings and precautions. As part of its continuous product development policy, Cordis reserves the right to change product specifications without prior notification. Please contact your Cordis representative for additional product availability information.

If you do not wish to receive any further emails from Cordis about our educational events and products/services, reply to this email with the word 'unsubscribe' in the subject line. More information about how Cordis processes personal data can be found at [https://www.cardinalhealth.co.uk/en\\_gb/search.html?q=privacy%20policy](https://www.cardinalhealth.co.uk/en_gb/search.html?q=privacy%20policy).

CORDIS, CORDIS Logo and MYNX CONTROL are trademarks of Cardinal Health and may be registered in the US and/or in other countries. © 2021 Cardinal Health. All Rights Reserved. 100581571 04/2021

