

Maj 2021

Hyvä asiakas,

Vaikka emme voi tehdä asiakaskäyntejä, haluamme pitää teidät ajan tasalla ja muistuttaa, että olemme tavoitettavissanne.

Lähetämme tämän uutiskirjeen mukana teille päivitettyä informaatiota verisuonten sulkijalaitteestamme (VCD), jota edustaa uuden sukupolven MYNX CONTROL<sup>TM</sup>. Laitteessa on muun muassa visuaalinen kireyden näyttö, mikä helpottaa toimenpidettä merkittävästi. Alempaa löydätte lisätietoja laitteesta.

Toivomme tuotteen kiinnostavan teitä. Annan mielelläni lisätietoja laitteesta ja vastaan kysymyksiinne. Mukavaa kevään jatkoa ja ollaan yhteydessä!

Ystävällisin terveisin,

Mikko Takkunen

Account Manager, Cordis

+358 40 519 77 96 / mikko.takkunen@cardinalhealth.com

## MYNX CONTROL<sup>TM</sup> Vascular Closure Device CLOSES WITH SECURITY. LEAVES WITHOUT A TRACE.



**MYNX CONTROL<sup>TM</sup>**  
VASCULAR CLOSURE DEVICE

**Closes with Security.  
Leaves without a Trace.**

**SECURE CLOSURE**   **SAFETY AND PATIENT COMFORT**   **EASE OF USE**

**Next-generation delivery system**  
• Deploy and retract: 2-button operation streamlines procedure  
• Sheath catch supports consistent deployment  
• Sealant is resorbed by the body within 30 days

**Tension indicator provides visual confirmation of device positioning for effective deployment**  
• Visual confirmation on handle  
• Correct balloon-to-wall apposition, for effective deployment

**Deploys through standard procedural sheath\***  
• No sheath exchange or up sizing required

Further information on MYNX CONTROL<sup>TM</sup> VCD [click here](#)

**Step 1: Position the Balloon**  
Result: Temporary Hemostasis

**Step 2: Deploy the Sealant**  
Result: Delivery

**Step 3: Remove the Device**  
Result: Secure Extravascular Closure

\*Note: MYNX CONTROL<sup>TM</sup> VCD is compatible with Medtronic Input<sup>®</sup> Introducer (11cm) sheaths, Cook Check-Flo<sup>®</sup> Performer<sup>®</sup> Introducer sheaths, and procedural sheaths longer than 12cm in effective length.

For Healthcare Professionals Only: Important information: Please refer to the instructions for use supplied with this device for indications, contraindications, side effects, suggested procedure, warnings and precautions. As part of its continuous product development policy, Cordis reserves the right to change product specifications without prior notification. Please contact your Cordis representative for additional product availability information.

© 2021 Cardinal Health and Cook. All other marks are the property of their respective owners.

CE 2797

Please Note: MYNX CONTROL<sup>TM</sup> VCD is incompatible with Medtronic Input<sup>®</sup> Introducer (11cm) sheaths, Cook Check-Flo<sup>®</sup> Performer<sup>®</sup> Introducer sheaths, and procedural sheaths longer than 12cm in effective length.

For Healthcare Professionals Only.

Important information: Prior to use, refer to the instructions for use supplied with this device for indications, contraindications, side effects, suggested procedure, warnings and precautions. As part of its continuous product development policy, Cordis reserves the right to change product specifications without prior notification. Please contact your Cordis representative for additional product availability information.

If you do not wish to receive any further emails from Cordis about our educational events and products/services, reply to this email with the word 'unsubscribe' in the subject line. More information about how Cordis processes personal data can be found at [https://www.cardinalhealth.co.uk/en\\_gb/search.html?q=privacy%20policy](https://www.cardinalhealth.co.uk/en_gb/search.html?q=privacy%20policy).

CORDIS, CORDIS Logo and MYNX CONTROL are trademarks of Cardinal Health and may be registered in the US and/or in other countries. © 2021 Cardinal Health. All Rights Reserved. 100581571 04/2021