

Hei.

Håper alt står bra til. Ettersom vi nå ikke kan besøke dere pga. den rådende situasjonen, er det viktig for oss å forsøke og holde dere oppdaterte.

Oversender dere her Cordis Nyhetsbrev om Mynx Control, vår nye generasjons VCD. En helt Extravaskulær lukkeanordning som er sikker og enkel å bruke. Håper at dere synes dette høres interessant ut. Dersom dere ønsker mer informasjon så ta kontakt med oss.

Ønsker dere en fortsatt fin vår.

Beste hilsener,
Monica Olsen
Account Manager, Cordis
+47 92 09 93 53 / monica.olsen@cardinalhealth.com

MYNX CONTROL™ Vascular Closure Device CLOSES WITH SECURITY. LEAVES WITHOUT A TRACE.



The brochure for MYNX CONTROL™ Vascular Closure Device features the Cordis logo and the tagline 'CLOSES WITH SECURITY. LEAVES WITHOUT A TRACE.' It highlights three key benefits: 'SECURE CLOSURE', 'SAFETY AND PATIENT COMFORT', and 'EASE OF USE'. The device is described as a 'Next-generation delivery system' with a '2-button operation' and 'simultaneous procedure'. It also mentions a 'Tension indicator provides visual confirmation of device positioning for effective deployment' and that it 'Deploys through standard procedural sheath'. The brochure includes three steps: 'Step 1: Position the Balloon Result: Temporary Hemostasis', 'Step 2: Deploy the Sealant Result: Delivery', and 'Step 3: Remove the Device Result: Secure Extravascular Closure'. A 'Further information on MYNX CONTROL™ VCD' link is provided. The CE mark is also visible.

Please Note: MYNX CONTROL™ VCD is incompatible with Medtronic Input® Introducer (11cm) sheaths, Cook Check-Flo® Performer® Introducer sheaths, and procedural sheaths longer than 12cm in effective length.

For Healthcare Professionals Only.

Important information: Prior to use, refer to the instructions for use supplied with this device for indications, contraindications, side effects, suggested procedure, warnings and precautions. As part of its continuous product development policy, Cordis reserves the right to change product specifications without prior notification. Please contact your Cordis representative for additional product availability information.

If you do not wish to receive any further emails from Cordis about our educational events and products/services, reply to this email with the word 'unsubscribe' in the subject line. More information about how Cordis processes personal data can be found at https://www.cardinalhealth.co.uk/en_gb/search.html?q=privacy%20policy.

CORDIS, CORDIS Logo and MYNX CONTROL are trademarks of Cardinal Health and may be registered in the US and/or in other countries. © 2021 Cardinal Health. All Rights Reserved. 100581571 04/2021