


Hej!

Hoppas att allt är bra. Eftersom vi inte kan besöka er pga. rådande situation så vill vi ändå försöka att hålla er uppdaterade. Översänder Cordis Nyhetsbrev med uppdaterad information om MYNX CONTROL™, vår nya generations VCD, en helt extravaskulär förslutningsanordning som är säker och enkel att använda.

Hoppas att ni tycker detta verkar intressant. Önskar ni mer information så hör gärna av er!  
Önskar er en trevlig fortsättning på våren!

Med bästa hälsningar,  
Christina Carlsson  
Account Manager Sweden South, Cordis  
+46 702 20 63 94 / christina.carlsson@cardinalhealth.com

## MYNX CONTROL™ Vascular Closure Device CLOSES WITH SECURITY. LEAVES WITHOUT A TRACE.



The brochure for MYNX CONTROL™ Vascular Closure Device features the Cordis logo and the tagline 'CLOSES WITH SECURITY. LEAVES WITHOUT A TRACE.' It highlights three key benefits: 'SECURE CLOSURE', 'SAFETY AND PATIENT COMFORT', and 'EASE OF USE'. The device is described as a 'Next-generation delivery system' with a 'simultaneous procedure' that includes a 'sheath catch' and a 'sealant' resorbed within 30 days. A 'Tension indicator' provides visual confirmation of device positioning for effective deployment. The device 'Deploys through standard procedural sheath' without the need for sheath exchange or up-sizing. The brochure includes three steps: 1. Position the Balloon (Result: Temporary Hemostasis), 2. Deploy the Sealant (Result: Delivery), and 3. Remove the Device (Result: Secure Extravascular Closure). A 'Further information on MYNX CONTROL™ VCD' link is provided. The CE mark is also visible.

Please Note: MYNX CONTROL™ VCD is incompatible with Medtronic Input® Introducer (11cm) sheaths, Cook Check-Flo® Performer® Introducer sheaths, and procedural sheaths longer than 12cm in effective length.

For Healthcare Professionals Only.

Important information: Prior to use, refer to the instructions for use supplied with this device for indications, contraindications, side effects, suggested procedure, warnings and precautions. As part of its continuous product development policy, Cordis reserves the right to change product specifications without prior notification. Please contact your Cordis representative for additional product availability information.

If you do not wish to receive any further emails from Cordis about our educational events and products/services, reply to this email with the word 'unsubscribe' in the subject line. More information about how Cordis processes personal data can be found at [https://www.cardinalhealth.co.uk/en\\_gb/search.html?q=privacy%20policy](https://www.cardinalhealth.co.uk/en_gb/search.html?q=privacy%20policy).

CORDIS, CORDIS Logo and MYNX CONTROL are trademarks of Cardinal Health and may be registered in the US and/or in other countries. © 2021 Cardinal Health. All Rights Reserved. 100581571 04/2021